KEIRSURGICAL

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KS 90 KEIRSURGICAL
Since 1923  Depuis 1923

Notre histoire

A family of innovative health care companies

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
STERILITY MAINTENANCE STUDY: DYNAMIC EVALUATION OF STERILIZED RIGID CONTAINERS AND WRAPPED INSTRUMENT TRAYS TO PREVENT BACTERIAL INGRESS.

PUBLISHED IN THE AMERICAN JOURNAL OF INFECTION CONTROL, DECEMBER 2015

RESEARCH CONDUCTED IN 2006 FOUND THAT STERILIZATION WRAP WAS MORE EFFECTIVE THAN RIGID CONTAINERS AT MAINTAINING STERILITY
SECOND STUDY DONE IN 10 YEARS TO HAVE SIMILAR FINDINGS²

HOW ARE CONTAINERS AND STERILIZATION WRAP USED?

Containers and sterilization wrap (also known as sterile packaging systems) protect surgical instruments from contamination from the time of sterilization until the surgical procedure, decreasing risk for SSIs.

THE METHODOLOGY¹

After being sterilized ...

... were tested for bacteria

THE TESTING PROTOCOL¹

Membranes in sterilization packaging

Sterilization occurs (Pre-Vac)

Aerosol chamber to emulate real-life environmental changes

Containers & wrapped trays wiped down


Membranes put into petri dishes & bacteria counted

THE RESULTS¹

Even unused rigid containers had high levels of contamination**

****72%** of unused containers showed bacterial ingress
*Out of 111 rigid containers tested, 14 (12.6%) had no bacterial ingress, 25 (22.5%) had ingress of 1-9 CFU, 52 (46.8%) had ingress of 10-99 CFU, and 20 (18.0%) had ingress >100 CFU.


100% of the wrapped trays maintained sterility until use

Research was funded by Halyard Health, Inc. (formerly Kimberly-Clark Health Care) and conducted by Applied Research Associates, an international research laboratory.

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Inspection Caisson– Surface externe

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≥ 1 leak: **test +**

no leak: **test -**

(1) FDS 98-053 : Test protocol to be performed on reusable sterile barrier systems (containers) to assess the watertightness of the container/lid closure (2014)

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
	Each Use	Yearly
Seals		
<ul style="list-style-type: none"> • Worn • Abraded • Cuts. Creases • Slipped from position • Lose • Hardened • Degrading (oxidized) • Stretched 	✓	
Dents, cuts, cracks	✓	
Handle	✓	
Filter Housing	✓	
Latches	✓	
Bolts, rivets, nuts for movement	✓	
Test bacterial barrier effectiveness		✓

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
Sommaire

Than KS

Keir Surgical Ltd.
126-408 East Kent Avenue South
Vancouver, BC V5X 2X7

T. 604.261.9596 / 800.663.4525 | F.
604.261.9549
info@keirsurgical.com
www.keirsurgical.com

ISO 13485:2003 Certified